



Make



every



day



amazing

Beltone **Amaze**™

Made for
iPhone | iPad | iPod

Works with
android

 **Beltone**™

80
YEARS
1940
2020

Beltone Amaze

The most complete hearing care solution now including a full family of hearing aids

Every day is filled with special moments that make our lives truly amazing. Every word, giggle or sound counts for something. Beltone Amaze ensures you won't miss a single one of those amazing moments.

Prepare to be amazed

The most complete and finest product in the industry, Beltone Amaze will give you an amazing hearing experience in all situations with a clearer, fuller and richer sound. Easy to use and offering full Beltone Remote Care whenever you need it, Beltone Amaze simply cannot be matched. Tap into excellent connectivity for crystal-clear streaming and personalized control and take advantage of the rechargeable battery* that will keep you going for a full day on a single charge.

Welcome to the world of Beltone Amaze.



*Only Beltone amaze RIE 63



A HEARING HOME RUN



See the Beltone Amaze video at
beltone.com/amaze



Amazingly natural sound

Enjoy a clearer, fuller and richer
sound experience

With Beltone Amaze, built on our latest and fastest micro-processing technology, you will experience a clearer, fuller and richer sound in any listening situation. Now you can hear all of life's beautiful and delicate sounds, from a grandchild's exciting story, a quick joke in a busy restaurant to birds singing in the distance. Get ready to experience more action, thrills or adventures when you are watching a movie. And, if you are a music or theatre lover, prepare to be amazed by the textured tones that make live performances deeper and more fulfilling.

Immerse yourself in layers of sound

Our industry-leading CrossLink Directionality 2 with Personal Sound ID offers unparalleled situational awareness. Conversations and ambient noise can now sound natural and you can select which sounds to focus on. And, if you get caught by a sudden loud noise, Beltone Amaze will adjust it to a more comfortable level.


Enjoy lush, more textured music

Listening to music puts a hearing aid to the ultimate test. With Beltone Amaze, the music has never sounded so rich and clear. Whether you are streaming your favorite melody or hearing it live at a symphony hall, Beltone Amaze will let you lose yourself in music – every time.

NEXT LEVEL SOUND QUALITY



See the Beltone Amaze sound
quality video at
beltone.com/amaze



**In a recent test, 95% of
listeners preferred Beltone
Amaze for music listening over
competitors' hearing aids***

*Jespersen et al, 2018

Amazingly **connected**

Excellent streaming quality
and easy control

Streaming music or phone calls directly to Beltone Amaze has never been more enjoyable. Experience the full action of your favorite movie streamed directly to your Beltone Amaze hearing aids, with no sounds missed.


Empowering personal control

You can easily control your hearing aids remotely with the **Beltone HearMax™ app**. With a few simple taps on the app, you can discreetly adjust your hearing aids directly from your iPhone®, iPad®, iPod touch®, Apple Watch® or an Android device.


From checking your battery status and adjusting programs and volume to using more advanced fine-tuning features that control speech focus, noise and wind filters, you will feel empowered and more satisfied with your listening experience.

A joy for you and for others

Imagine enjoying a great movie on TV with your family – but at your own volume. Or having a conversation in the car without missing a word. With the help of Beltone's wireless accessories you will be able to catch every conversation and watch TV with your loved ones at your own volume.



When comparing sound quality to competitor hearing aids for TV streaming, listeners preferred Beltone Amaze 71% of the time



When comparing to competitor hearing aids that stream from iPhone, listeners chose the sound quality of Beltone Amaze 64% of the time*



*Jespersen et al, 2018

Amazingly personalized care

With you every step of the way


Getting used to your hearing aids can take time, and your hearing needs can change over time. Beltone Remote Care enables easier communication between you and your hearing care professional, so your hearing aids can be continuously adjusted to your specific needs. This makes your journey toward better hearing smoother and easier, every step of the way.

Easy, convenient care

With the help of the free Beltone HearMax app, you can easily send a request to your hearing care professional and receive the new settings with a simple tap on your phone, without the need for an office visit. Beltone Remote Care will save you time and give you confidence that help is just a tap away.

Beltone Remote Care

is easy as 1-2-3




1 You request assistance via the HearMax app.



2 Your hearing care professional reviews and sends the new settings back to your app.



3 Tap on your app screen to accept the update. You can also rate the new settings in the app.



Beltone Remote
Care App



PROVIDING CARE. ANYWHERE.



See the Beltone Amaze
Remote Care video at
beltone.com/amaze

Amazingly rechargeable

The freedom to charge into life

Make battery changes a thing of the past

Thanks to Beltone Amaze's integrated lithium-ion battery, you no longer need to replace batteries*. And with the most user-friendly portable charger in the industry, you can charge confidently into every day.


Up to 30 hours on one charge

Beltone Amaze will last you a long time, up to 30 hours, when fully charged. Recharging is simple - place your hearing aids in the charger overnight and they are ready again in the morning. If you forgot to charge your hearing aids, no worries – in just 30 minutes of charging Beltone Amaze will be ready to go another 8 hours before the next charge. Extremely durable – even after 4 years of use – the battery in your Beltone Amaze will retain 80% of its capacity.



Portable charger and carrying case

Your Beltone Amaze hearing aids comes with a sleek portable charger that also serves as a carrying case. You will find using the charger to be amazingly simple – the LED lights will show you the current battery level and if your hearing aids were inserted correctly and charging. The charger itself has an internal battery with an impressive capacity to give you 3 full charges.

Battery life 	Charging required
30 hours	3 hours
16 hours	1 hour
8 hours	30 minutes

Beltone Amaze

Product family

CIC



Completely-in-the-Canal

The first CIC with 2,4 GHz streaming. The smallest and most discrete of the Amaze custom instruments. Designed for individuals with mild to moderate hearing loss.

ITC



In-the-Canal

The ITC is the medium size custom option where the instrument sits slightly out of the ear canal and into the outer ear. This size allows more options to be added. Designed for individuals with mild to moderately severe hearing loss.

ITE



In-the-Ear

The largest of the custom instruments which fills the outer ear portion of the ear. The ITE is the most flexible size for adding options. The larger size also allows more powerful components inside the hearing aid to help patients with mild to severe hearing loss.

MIH



Microphone-in-Helix

This style is exclusive to Beltone. The microphone is taken out of the body of the instrument and placed in the fold of the ear above the canal. The placement of the microphone helps block wind noise and enhance sound quality. Designed for active individuals with mild to moderately severe hearing loss.

RIE



Receiver-in-the-Ear

One of the most popular styles of hearing aids. The hearing aid rests over the ear, but the speaker (receiver) is placed in the ear canal via a dome or custom ear mold. The receiver is connected to the hearing aid by a thin wire that is very discrete. The placement of the receiver in the canal provides acoustic and sound quality benefits. Designed for individuals with mild to severe hearing loss.

BTE



Behind-the-Ear

Easy-to-wear BTE devices are small, and feature a curved casing designed to nestle comfortably behind the ear. By matching flesh tone or hair color, they hide with ease. A transparent tube connects the device to a custom-created ear mold worn in your ear canal. Ear molds are lightweight and comfortable. Designed for individuals with mild to severe hearing loss.

* only available on select models

Beltone Amaze

2020 BIG Innovation Award Winner



Breakthrough technology allowing direct streaming to Android operating system devices and preserving battery life via Bluetooth Low Energy.

© 2020 GN Hearing Care Corporation. All rights reserved. Beltone is a trademark of GN Hearing Care Corporation. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google LLC. The Android robot is reproduced or modified from work created and shared by Google and used according to terms described in the Creative Commons 3.0 Attribution License. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.

Manufacturer According
to FDA:
Beltone
8001 E Bloomington Fwy.
Bloomington, MN 55420
U.S.A.
(800) 621-1275
beltone.com