

Premier Hearing Care

Hear. There. Anywhere.



Beltone Trust™



Made for
 iPod  iPhone  iPad

The world's most complete **hearing care** experience

Life is full of pleasant surprises! From having fun with friends and family to travelling the world for business or pleasure, it should be enjoyed. And your hearing should never cause you to miss a single experience.

Beltone Trust prepares you for life's special moments better than any other hearing aid, wherever you are. It offers exceptional sound, discreet and customized design and personal control. You deserve the best, and Beltone Trust brings the premier hearing care you need to a whole new level. The groundbreaking new Remote Care™ feature gives you extra support in-between scheduled appointments with your hearing care professional, no matter where you are.

Beltone Trust is your trusted partner that that helps you get the most out of life – hear, there or anywhere you are.

WATCH THE MOVIE

Learn how Beltone Trust can change your quality of life.



See the Beltone Trust movie at beltone-hearing.com



WATCH THE MOVIE

Learn how Personal Sound ID and CrossLink Directionality 2 help you hear the way you were born to.



See the Beltone Trust movie at beltone-hearing.com

Trust in the most **personalized** and natural hearing

Our ears work together in a wonderful way known as binaural hearing. This is when our brains receive information from each of our ears, and empower us to listen selectively and understand our surroundings while hearing. Average hearing aids find it difficult to imitate this, but Beltone Trust includes two groundbreaking innovations that help you hear in this more natural way: Personal Sound ID™ and Crosslink Directionality™ 2.

More natural hearing with Personal Sound ID

This unique feature helps restore the natural sound level differences between your ears. You benefit from a richer sound palette and the ability to hear where sounds originate from.

A complete sound experience with CrossLink Directionality 2

CrossLink Directionality 2 offers you truly individual listening by monitoring your sound environment and automatically switching to the most natural speech focus settings. It combines with Personal Sound ID to help you to better recognize voices and the direction sounds come from. And you'll understand that speech better, too. According to a recent study*, Beltone hearing aids make it significantly easier to understand speech than major competitors. This is especially true when you want to hear someone talking behind you or to your side.

*Source: Study by GN Group, Jespersen and Kirkwood, 2016

Trust in great sound

Your Beltone Trust hearing aids adapt to you and your lifestyle. You'll enjoy a new level of hearing quality in a wider range of situations.

Smart Gain Pro™ adjusts volume, adapting automatically to changes in your environment. Whether you are in a noisy, crowded restaurant or a quiet park, your hearing will be at its best.

Feedback Eraser™ with Whistlestop™ cancels feedback for your comfort. This means you don't have to endure unpleasant noises and whistling that can sometimes occur when you use your phone or hug a loved one or friend.

Ear to Ear synchronization means that you can easily move from a quiet place to a noisier one. This feature monitors your environment and automatically coordinates your hearing aids, so you always experience a smooth sound transition.



Trust in personal control

The **Beltone HearMax** app gives you an easy and discreet option to control your hearing aids. Simply download the free app and start adjusting volume, listening programs and much more directly from your iPhone, iPad, iPod touch, Apple Watch or the most popular Android phone models. Access it all from an new, intuitive one-tap main screen.



The Beltone HearMax app lets you change volume, programs and much more directly from your smartphone or tablets.



WATCH THE MOVIE

See how the Beltone HearMax app lets you control your hearing aids. It is as easy as sending a text message or email!



See the Beltone Trust movie at beltone-hearing.com



Trust that you'll always stay connected

Direct audio streaming (Made for iPhone)

Innovative wireless technology enables you to stream sound directly from your iPhone, iPad or iPod touch to eight different Beltone Trust hearing aid models, without having to wear any sort of adapter around your neck. So you can hear voices from phone calls, listen to your favorite music or follow the voice directions of a GPS – all through your Beltone Trust hearing aids.

WATCH THE MOVIE

See how you can stream sound directly from an iPhone, iPad or iPod touch to your hearing aids. It's a simple connection that acts like wireless headphones.




See the Beltone Trust movie at beltone-hearing.com



Trust in premier service


Beltone Remote Care gives you the most personalized hearing care experience possible. This exclusive, revolutionary technology allows your hearing aids to be fine-tuned wherever you are. You simply send a request via the Beltone HearMax app and your hearing care professional will send new settings you can download with your mobile phone.

Beltone Remote Care makes it possible for you to receive hearing aid adjustments from your hearing care professional directly via your iPhone or Android phone – wherever you are.



WATCH THE MOVIE

See how your hearing care professional can use Remote Care technology to fine-tune your hearing aids anywhere, giving you the ultimate in hearing confidence – wherever you are.



See the Beltone Trust movie at beltone-hearing.com

Trust in discreet design

With Beltone Trust, you'll always look your best and feel comfortable. You can choose from a variety of Custom, Receiver-in-the-ear and Behind-the-ear models that come in an assortment of colors.

Beltone Trust is small, virtually invisible and built to last. The HPF⁸⁰ NanoBlock coating blocks moisture and dirt and the Wind Noise Reduction feature makes it easy to enjoy breezy days. So if whether you're chatting with friends indoors or outside in tough weather, you don't need to worry about your hearing aids.



*Invisible-in-the canal
(IIC) custom model*



*Receiver-in-the-ear
(RIE) 63*



IIC, CIC, ITC, ITE, MIH-S, MIH, 63 RIE, 64 RIE, MINI BTE 66, 76 BTE, 86 BTE



Beltone Trust Receiver-in-the-Ear (RIE) and Behind-the-ear (BTE) colors

Beltone. Helping the world hear better

Beltone was founded on the act of helping a friend to enjoy life more. Since 1940, we have provided knowledge, tools, service and training to the professionals servicing the hearing impaired. Beltone makes quality hearing care available to more people. To do so, we develop technically optimal hearing solutions without ever losing focus on the individual's needs. It is our belief that no individual should be denied the possibility of an improved life with better hearing.

© 2017 GN Hearing Care Corporation. All rights reserved. Beltone is a trademark of GN Hearing Care Corporation. Apple, the Apple logo, iPhone, iPad, iPod touch and Apple Watch are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc., registered in the U.S. and other countries. Android, Google Play and the Google Play logo are trademarks of Google Inc.

For product compatibility, please check
www.beltone-hearing.com/compatibility

Worldwide Headquarters

Beltone A/S
Lautrupbjerg 7
DK-2750 Ballerup
Denmark
Tel.: + 45 45 75 11 11
Fax: + 45 45 75 11 19
Web: www.beltone-hearing.com

