

### **Beltone. Helping the world hear better**

Beltone was founded on the act of helping a friend to enjoy life more. Since 1940, we have provided knowledge, tools, service and training to the professionals servicing the hearing impaired. Beltone makes quality hearing care available to more people. To do so, we develop technically optimal hearing solutions without ever losing focus on the individual's needs. It is our belief that no individual should be denied the possibility of an improved life with better hearing.

Beltone is a trademark of GN Hearing Care Corporation. © 2015 GN Hearing Care Corporation. All rights reserved.

\*SmartRemote app is compatible with iPhone running iOS 7.X or later. Apple, the Apple logo and iPhone are trademarks of Apple Inc., registered in the U.S. and other countries. App Store is a service mark of Apple Inc. The Android version of the Beltone SmartRemote app is available on select Android devices running version 2.3.5 or later. Android and Google Play are trademarks of Google Inc. Dolby and the double-D symbol are registered trademarks of Dolby Laboratories.

#### **Worldwide Headquarters**

Beltone A/S  
Lautrupbjerg 7  
DK-2750 Ballerup  
Denmark  
Tel.: + 45 45 75 11 11  
Fax: + 45 45 75 11 19  
Web: [www.beltone-hearing.com](http://www.beltone-hearing.com)



# Beltone **Direct Line™**

## Extend your hearing with **wireless accessories**



# Experience true hearing freedom

Beltone Direct Line™ wireless accessories make it easier for you to engage and interact in any hearing situation and link to the world around you. Whether you are watching television, following a conversation at a distance or talking on the phone, we have a wireless accessory that can help you hear better.

Beltone Direct Line features Beltone's own 2.4 GHz wireless system. As well as providing you with a better quality audio signal, this proven technology allows you to stream stereo sound directly into your hearing aids – without bothersome connectors around your neck. All you need are a few handy devices and you're right back in the conversation, enjoying TV or talking on the phone.

## Direct TV Link 2

Enjoy streaming high-quality stereo sound from your TV, stereo and PC directly to your hearing aids.

## Direct Remote Control 2

Conveniently and discretely adjust your hearing aids and control the audio volume streamed from your Direct TV Link 2, Phone Link 2 and myPAL.

## Direct Phone Link 2

Experience crystal clear phone conversations and enjoy the freedom of basic remote control functionality.

## myPAL Micro (Personal Audio Link)

Hear clearly at a distance as myPAL Micro streams almost any sound you like directly into your hearing aids. Great for noisy places.

## myPAL Pro (Personal Audio Link)

Besides the features of the myPAL Micro, the myPAL Pro provides you with a clever table mic functionality which automatically enables it to pick up the voices of multiple speakers around the unit when you place it horizontally. Perfect for dinners with friends and family or a business meeting.



# Simply extend your hearing range

Imagine the feeling of being right next to the person you want to hear – even though you're over 20 meters away – that's the kind of hearing experience the new Beltone myPAL (Personal Audio Link) offers. Clip it on a tour guide in a museum, your grandchild in the car or a good friend in a cafe – and get close to what you want to hear.

Non face-to-face communication can be a challenge for everyone – hearing aids or not. The Beltone myPAL is a small portable microphone you hand over to the person you want to hear. It captures the sound clearly, even in noisy environments, and streams it directly into your hearing aids, bringing the conversation or entertainment closer to you. It's that simple.

## Beltone myPal Micro

- **Wireless range of up to 25 m (clear line of sight) keeps you in touch even at a distance**
- **Directional microphone for clear, crisp sound**
- **On-board volume control and mute functionality**
- **Drop detection, which mutes the sound transmission if it's dropped from a height above 75 cm (29.5 inches)**
- **Visual status indicators help you keep an eye on current status and mode of operation**
- **Battery life of 10 hours for all-day confidence**

## Beltone myPAL Pro

*All the features of myPAL Micro, plus:*

- **When used horizontally myPAL Pro automatically switches into a table mode optimized for picking up the voices of multiple speakers**
- **Extended choice of connectivity options**



**NEW**



## The **ultimate TV** experience for everyone

Remember when you turned up the volume to hear the TV or stereo? Now you and your family can sit together in the room, with the TV volume at a comfortable level for everyone and you as a fully engaged conversation partner.

With the Beltone Direct Line TV Link 2, you can enjoy clear sound from your TV or any other audio device – at your own preferred volume – and still participate in conversations around you. It connects directly to virtually any audio device: TVs, stereos or computers, and delivers entertainment directly to your hearing aids.



### Direct TV Link 2

- **Streams stereo sound from TV or any other audio source to both hearing aids**
- **Streams sound at least 7 meters with clear line of sight**
- **Reconnects automatically if you get out of range and come back within 5 minutes**
- **Works with both digital and analogue audio inputs**
- **Accepts Dolby® Digital encoded audio inputs**

# Full control at your fingertips

Controlling the way you choose to hear has never been easier. The intuitive Beltone Direct Line Remote Control 2 makes hearing simple. You can adjust the volume in one or both hearing aids, change programs, or easily switch between sound streamed from your TV, radio, myPAL or PC. The large LCD display shows your programs, volume and battery level at a glance.

Now you can even adjust the balance between the sounds around you and the streaming audio to make sure you hear everything you want.

## Direct Remote Control 2

- **Adjusts volume on one or both hearing aids**
- **Changes hearing aid programs**
- **Switches between streaming from your TV, radio or laptop**
- **Balances the hearing aid volume with the volume of your streaming devices**
- **Blocks out surrounding sounds with the Mute function**
- **Offers discreet and clear overview of your hearing aids settings**



# Clearer phone conversations

Hear the phone ring directly in your hearing aids, even when your phone is put away in a pocket or a handbag. Speak hands-free, and enjoy clear conversation in both ears. Use the mute button to eliminate background noise when you're talking in loud surroundings.

The Phone Link 2 also serves as a remote control to change your environmental programs and adjust hearing aid volume whether or not you're on the phone. And the Beltone SmartRemote™ app works with the Phone Link 2 to turn an iPhone® or Android™ phone into a discreet, full-functioning remote control.\*

## Direct Phone Link 2

- **Strong, clear speech signal**
- **Mute and Program buttons offer a convenient way to control basic hearing aid functionalities**
- **Dual microphones and a noise reduction feature help make speech easier to understand**
- **Stream high quality stereo audio from your phone**



SmartRemote App



Phone Link 2

